

**Intended Use**

This reagent is intended for *in vitro* quantitative determination of SGPT in serum or plasma.


- IFCC recommended methodology
- Linear up to 1000 U/L
- Working reagent can be prepared as per requirements

**Clinical Significance**

It is present in most of the tissues, but mainly found in the liver. Increased levels are found in hepatitis, cirrhosis, obstructive jaundice & other hepatic disease. SGPT activity is markedly elevated even before clinical signs of jaundice become apparent in disease associated with hepatic necrosis. Slight elevations are also found in myocardial infarction.

**Principle**

Kinetic determination of Alanine Aminotransferase (ALT) according to the following reaction.

**Kit Components**

Reagent/Component	Product code 11409005	Product code 11409003	Product code 11409006	Description
SGPT (S.L) R1	2 x 24 mL	3 x 40 mL	4 x 100 mL	Tris buffer (pH 7.5) 110 mmol/L L-Alanine 600 mmol/L LDH > 1500 U/L
SGPT (S.L) R2	2 x 6 mL	3 x 10 mL	4 x 25 mL	alpha -ketoglutarate 16 mmol/L NADH 0.24 mmol/L

**Risk & safety**

Material Safety data sheets (MSDS) will be provided on request

**Reagent Preparation**

Mix 4 volume of Reagent 1 (R1) with 1 volume of Reagent 2 (R2)

This Working reagent is stable for 30 days at 2-8°C.

Note: Discard the working reagent, if the blank absorbance is less than 1.0 at 340 nm.

**Reagent Storage and Stability**

The sealed reagents are stable up to the expiry date stated on the label, when stored at 2- 8°C and protected from light.

**Open Vial Stability**

Once opened, the reagent is stable up to 4 weeks at 2-8°C if contamination is avoided.

**Onboard Calibration Stability**

On-board Calibration stability is 20 days

**Reagent Deterioration**

Turbidity or precipitation in any kit component indicates deterioration and the component must be discarded. Values outside the recommended acceptable range for the Agappe Qualicheck Norm & Path control may also be an indication of reagent instability and associated results are invalid. Sample should be retested using a fresh vial of reagent.

**Precaution**

To avoid contamination, use clean laboratory wares. Use clean, dry disposable pipette tips for dispensing. Close reagent bottles immediately after use. Avoid direct exposure of reagent to light. Do not blow into the reagent bottles.

This reagent is only for IVD use and follow the normal precautions required for handling all laboratory reagents

**Waste Management**

Reagents must be disposed off in accordance with local regulations

**Sample**

Serum /Plasma (free of haemolysis)

**Interferences**

No interference for

- Bilirubin up to : 10 mg/dL
- Ascorbic acid up to : 500 mg/dL
- Haemoglobin up to : 1000 mg/dL

**Materials provided**

SGPT Reagent R1 & R2

**Materials required but not provided**

- Pipettes & Tips
- Test Tubes & racks
- Timer
- Incubator
- Analyzer

**Test Parameters**

	Normal procedure	High Linearity procedure
Mode of Reaction	Kinetic	Kinetic
Slope of reaction	Decreasing	Decreasing
Wavelength	340 nm	340 nm
Temperature	37°C	37°C
Factor	1745	1745
Linearity	350 U/L	1000 U/L
Blank	DI Water	DI Water
Delay	60 sec	60 sec
No of reading	3	3
Interval	60 sec	20 sec
Sample volume	100 µL	100 µL
Reagent volume	1000 µL	1000 µL
Cuvette	1 cm light path	1 cm light path

Application parameters for various instrument are available. Please contact customer support department for specific information.

**Unit Conversion**

Traditional Unit	SI Unit	Conversion from Traditional to SI
U/L	µKat/L	x 0.017

**Calibration**

Agappe multicalibrator is recommended for calibration of this assay on fully auto analyzer.

Use provided factor (1745) for estimation of this assay on semi auto analyzer

**Procedure notes**

Laboratory procedure for Semi Auto Analyzer	
Working reagent	1000 µL
Sample	100 µL
Mix and incubate at 37°C for 1 minute. Measure the change in absorbance per minute ( $\Delta$ OD/min) during 3 minutes.	
High Linearity Procedure	
Mix and incubate for 1 minute at 37°C. Read the change in absorbance per 20 sec during 1 minutes.	

**Calculation**

SGPT activity (U/L) = ( $\Delta$  OD/min) x 1745

**SYMBOLS USED ON THE LABELS**

IVD IN VITRO DIAGNOSTIC USE SEE PACKAGE INSERT FOR PROCEDURE LOT NUMBER MANUFACTURER'S ADDRESS MANUFACTURING DATE EXPIRY DATE TEMPERATURE LIMIT

**AGAPPE DIAGNOSTICS LTD.**

'Agappe Hills', Dist. Ernakulam, Kerala, India-683 562.  
Tel. +91 484 2867 000 | Customer Support No.: 1800 425 7151 (Toll Free)  
customersupport@agappe.in | www.agappe.com

REV. NO.: ADL/IFU/SGPT/LIQ/R02

ISO 9001 : 2015  
EN ISO 13485: 2016

**Quality control**

It is recommended to use Agappe Qualicheck Norm & Path ( 11601001) to verify the performance of the assay. Each laboratory has to establish its own internal quality control scheme and procedure for corrective action, if control do not recover within the acceptable range.

**Reference Range**

It is recommended that each laboratory has to establish its own reference values.

The following value may be used as guide line.

Serum up to : 49 U/L

Results obtained for patient samples are to be correlated with clinical findings of patient for interpretation and diagnosis.

**Performance****1. Linearity**

This reagent is linear up to 1000 U/L.

If the concentration is greater than 350 U/L, follow the high linearity procedure to get higher linearity of 1000 U/L.

If the concentration is greater than linearity dilute the sample with normal saline and repeat the assay. Multiply the result with dilution factor.

**2. Comparison**

A comparison study has been performed between Agappe reagent and another internationally available reagent yielded a correlation coefficient of  $r^2= 0.9997$  and a regression equation of  $y = 0.9906x$ .

**3. Precision**

Intra Run		
	Control Level 1	Control Level 2
n	20	20
Mean (U/L)	38.0	115.9
SD	1.08	3.29
CV(%)	2.86	2.84

Inter Run		
	Control Level 1	Control Level 2
n	20.00	20.00
Mean (U/L)	38.7	115.1
SD	1.2	2.7
CV(%)	3.1	2.5

Accuracy (U/L)		
Control	Expected Value	Measured Value
Control Level 1	28 ± 10.2	31.3
Control Level 2	98.5 ± 25	106.1
Qualicheck Norm	48 ± 9.8	52.4
Qualicheck Path	125 ± 30	130

**4. Sensitivity**

Lower detection Limit is 0.5 U/L

**Bibliography**

- Clin. Chem, Acta. 105, 147-172 (1980)
- Thefeld, W., *et al.*; Dtsch. Med Wschr.99, 343 (1994)

**SYMBOLS USED ON THE LABELS**

**IVD** IN VITRO DIAGNOSTIC USE  SEE PACKAGE INSERT FOR PROCEDURE **LOT** LOT NUMBER  MANUFACTURER'S ADDRESS  MANUFACTURING DATE  EXPIRY DATE  TEMPERATURE LIMIT

**AGAPPE DIAGNOSTICS LTD.**

'Agappe Hills', Dist. Ernakulam, Kerala, India-683 562.  
Tel. +91 484 2867 000 | Customer Support No.: 1800 425 7151 (Toll Free)  
customersupport@agappe.in | www.agappe.com

REV. NO.: ADL/IFU/SGPT/LIQ/R02

**CE** ISO 9001 : 2015  
EN ISO 13485:2016