

Intended Use



This reagent is intended for *in vitro* quantitative determination of urea in serum, plasma & urine.

- Urease / GLDH methodology
- Linear up to 300 mg/dL

Clinical Significance

Proteins cannot be stored in human body, so excess should be broken down. Amino acids which form the components of proteins break down to give ammonia. This is toxic & so through a series of chemical reactions (urea cycle) non toxic urea is produced & this is released into the blood which is filtered in the kidney & excreted in the urine.

Elevated levels are seen during increased protein breakdown, dehydration, vomiting, and diarrhea. Also seen in any kind of renal disorder like Glomerular nephritis, chronic nephritis & nephritic syndrome.

Principle

Kit Components

Reagent/Component	Product Code 12006023	Description
Urea U V R1	4 x 50 mL	Buffer (pH 7.6) 100 mmol/L ADP 0.7 mmol/L α -ketoglutarate 9.0 mmol/L
Urea U V R2	2 x 30 mL	GLDH \geq 1100 U/L Urease \geq 6500 U/L NADH 0.25 mmol/L 2-Oxoglutarate 5 mmol/L

Risk & Safety

Material Safety data sheets (MSDS) will be provided on request

Reagent Preparation

Urea U V R1 & R2 Reagents are ready to use.

Reagent Storage and Stability

The sealed reagents are stable upto the expiry date stated on the label, when stored at 2-8°C.

Open Vial Stability

Once opened the reagents are stable up to 60 days if contamination is avoided.

On-board Calibration Stability

Calibration is stable for 20 days.

Reagent Deterioration

Turbidity or precipitation in any kit component indicates deterioration and the component must be discarded. Values outside the recommended acceptable range for the Agappe Qualicheck Norm & Path control may also be an indication of reagent instability and associated results are invalid. Sample should be retested using fresh vial of reagent.

Precaution

To avoid contamination, use clean laboratory wares. Close reagent bottles immediately after use. Avoid direct exposure of reagent to light. Do not blow into the reagent bottles.

This reagent is only for IVD use and follow the normal precautions required for handling all laboratory reagents.

Waste Management

Reagents must be disposed off in accordance with local regulations.

Sample

Serum, Plasma (Do not use lipemic or hemolysed sample) - Do not use anticoagulants containing fluorides and ammonium ions.

Urine - Diluted to 1/100 with distilled water. Multiply the result with dilution factor.

Interferences

No interference for
Bilirubin up to 20 mg/dL
Ascorbic acid up to 50 mg/dL
Haemoglobin up to 1000 mg/dL

Materials provided

Urea U V R1 & R2 Reagent

Reagents required but not provided

Multicalibrator (Product Code: 11610001), Qualicheck Norm (Product Code: 11601003), Qualicheck Path (Product Code: 11601002)

Unit Conversion

Traditional Unit	SI Unit	Conversion from Traditional to SI
mg/dL	mmol/L	x 0.1665

Calibration

Agappe Multicalibrator (Product Code: 11610001) is recommended for calibration of the assay.

Quality Control

It is recommended to use Qualicheck Norm (Product Code: 11601003) or Qualicheck Path (Product Code: 11601002) to verify the performance of the measurement procedure.

Each Laboratory has to establish its own internal quality control scheme and procedures for corrective action if controls do not recover within the acceptable tolerance.

Reference Range

It is recommended that each laboratory should establish its own reference values.

The following value may be used as guide line.

Serum/Plasma : 10-50 mg/dL

Urine : 20-35 gm/24 hr

Results obtained for patient samples are to be correlated with clinical findings of patient for interpretation and diagnosis.

Performance
1. Linearity

The reagent is linear up to 300 mg/dL. If the concentration is greater than linearity (300 mg/dL), dilute the sample with normal saline and repeat the assay. Multiply the result with dilution factor.

2. Comparison

A comparison study has been performed between Agappe reagent and another international available reagent yielded a correlation coefficient of $r^2 = 0.9828$ and a regression equation of $y = 0.9706x$.

3. Precision

Control	Intra Run		Inter Run	
	Level 1	Level 2	Level 1	Level 2
n	20	20	20	20
Mean (mg/dL)	36	96	35.8	96.1
SD	1.2	2.0	1.0	2.4
CV(%)	3.3	2.08	2.7	2.4

Accuracy (mg/dL)

Control	Expected Value	Measured Value
Control Level 1	32.6 \pm 9.8	32.3
Control Level 2	102 \pm 23	103
Qualicheck Norm	40 \pm 8.42	41.6
Qualicheck Path	126 \pm 18	127.3

4. Sensitivity

Lower detection Limit is 1 mg/dL

Bibliography

1. Talke, H. and Schubert G.E. Kiln-Wocchsr 43: 174 (1965)
2. Kassirer, J.P., New Eng. J. Med. 285, 385 (1971)

SYMBOLS USED ON THE LABELS

 IN VITRO DIAGNOSTIC USE
  SEE PACKAGE INSERT FOR PROCEDURE
  LOT NUMBER
  MANUFACTURER'S ADDRESS
  MANUFACTURING DATE
  EXPIRY DATE
  TEMPERATURE LIMIT

Accute 40 FR Assay Parameter

Page 1	
Test	UREA U.V
Reaction Mode	Rate
Reference Test ID	**
Test WL	340
Blank WL	NA
Test Read Timing	36-42
Blank Read Timing	NA
Sample	3 µL
R1	200 µL
R2	50 µL
Stirrer	ON
Cal Mode	std
Page 3	
Calibration Mode	Linear
** Not applicable	
# - Input Programme Number only for sample Blanking parameters	

SYMBOLS USED ON THE LABELS

IVD IN VITRO DIAGNOSTIC USE  SEE PACKAGE INSERT FOR PROCEDURE **LOT** LOT NUMBER  MANUFACTURER'S ADDRESS  MANUFACTURING DATE  EXPIRY DATE  TEMPERATURE LIMIT



AGAPPE DIAGNOSTICS LTD.

'Agappe Hills', Dist. Ernakulam, Kerala, India-683 562.
 Tel. +91 484 2867 000 | Customer Support No.: 1800 425 7151(Toll Free)
 customersupport@agappe.in | www.agappe.com

REV. NO.: ADL/IFU/UUV/40FR/R01

CE ISO 9001:2015
 EN ISO 13485:2016