



## BIOLOGICAL MICROSCOPE- INCLINED, MONOCULAR WITH LED

CAT. No. **LT- 327**

Achromatic Objective: 5x, 10x, 45x (SL),  
100x (SL) Oil immersion, Eye Piece 10x W.F.  
Quadruple revolving nose piece, Coaxial mechanical stage.  
Graduated precise slow motion.  
Reading to 0.002mm condenser unit N.A. 1.25 with  
Iris diaphragm with LED light.  
Complete with thermocol box + wooden box.



## BIOLOGICAL MICROSCOPE- BINOCULAR WITH LED

CAT. No. **LT- 328**

Achromatic Objective: 5x, 10x, 45x (SL), 100x (SL)  
Oil immersion, Eye Piece 10x W.F. (paired)  
Quadruple revolving nose piece., Coaxial mechanical stage.  
Binocular head with hard coated prisms.  
Graduated precise slow motion, with LED light.  
Complete with thermocol box + wooden box.

## MICROSCOPIC SPARE PARTS



## ACCESSORIES FOR MICROSCOPE

CAT. No.

**LT- 329**

Mechanical Stage, Condenser Unit with Iris Diaphragm, Objectives - 100x, 45x, 40x, 10x, 5x etc.  
Eye pieces WF - 10 x, 5x, 15x, Binocular head, Dual head attachment, Plano concave reflector 50mm etc.