

TUBERCULOSIS

TUBERCULIN PPD

For *in vivo* Detection of Sensitisation by *Mycobacterium tuberculosis*

- Diluted, ready to use Purified Protein Derivative solution for Mantoux test
- Highly purified material, calibrated against WHO approved standard preparation (RT 23) for better accuracy
- Available in 5 mL and 1 mL pack size
- Filled in amber coloured glass vials to prevent the loss of potency in regular usage and stabilized with Tween 80 to prevent absorption of tuberculin proteins to glass surface

FORTH COMING

Crystal[®] TB Confirm

Immuno-chromatographic Rapid Visual Test for Detection of growth of Tubercle bacilli in Culture

- Differential detection of *M. tuberculosis* complex (MTC) from Non-Tuberculous Mycobacteria (NTM)
- Detection is based on two separate targets MPT64 and MTCAg. Use of additional target - MTCAg, ensures detection of mutants not secreting MPT64 antigen

RHEUMATOLOGY

RA / CRP / ASO

Latex Agglutination Slide Test for Qualitative and Semi-Quantitative Detection of Rheumatoid Factor (RF) C-Reactive Protein (CRP) Anti Streptolysin-O (ASO)

- Rapid and simple test procedure
- Polystyrene latex particles of optimized size are used
- Serum dilution not required
- All the components & accessories required to perform test are included in the kit

ANTIBIOTIC DISC

CLAIRO[™] COMBI

Combination of Antibiotic Discs for Sensitivity Test

CLAIRO[™] MONO

Single Antibiotic Discs for Bacterial Sensitivity Test

TYPHOID

WIDAL ANTIGEN SET (Slide / Tube)

Stained Salmonella Antigens for Slide & Tube Tests

- Antigens are prepared from specially selected smooth, non-fimbriate strains
- Uncompromised sensitivity and accuracy with uniform clump formation
- Multiple pack sizes and long shelf life (18 months)
- Reagents are supplied in USP type 1 amber colored glass vial and offered in different combinations
 - 2 Antigens set - O & H
 - 4 Antigens set - O, H, AH & BH
 - 8 Antigens set - O, H, AH, AO, BO, BH, CO & CH
- Also available in Tube Test

CHOLERA

Crystal[®] Vc

Rapid Visual Immuno-chromatographic Test for Detection of *V. cholerae* in Stool

- Rapid, visual qualitative immunoassay with inbuilt control line
- Two separate bands for simultaneous & differential detection of *V. cholerae* O1 & O139 in a single test
- All the components & accessories required to perform test are included in the kit

