

# AGAPPE Mispa HX88

## 6 Part Hematology Analyzer with RET + IPF

---

### Why Choose?

The Agappe Mispa HX88 stands out for its advanced technology, robust build, and high diagnostic accuracy, making it an ideal solution for busy hospital and diagnostic laboratories. This analyzer brings together quality, speed, and operational simplicity for reliable healthcare delivery.

- Delivers fast, high-precision results with advanced laser flow cytometry.
  - Reduces manual handling—safeguarding both patient and technician.
  - Flexible operational modes for any laboratory need.
  - Processes up to 100 samples per hour for busy labs.
  - Real-time STAT handling for emergency diagnostics.
  - User-friendly touchscreen and smart interface ensure easy operation.
  - Powerful abnormal cell detection with 37+ parameters.
  - Minimal maintenance and great durability for round-the-clock use.
- 

### Unique Selling Points

- Third-generation nucleic acid fluorescence staining ensures ultra-precise cell differentiation.
  - Tri-angle laser scatter flow cytometry provides high accuracy for abnormal cell detection (including NRBC, IG, RET, IPF).
  - Up to 100 tests per hour to support high-volume workloads.
  - Multiple operational modes for versatile sample processing: whole blood, low-value leukocyte, predilution, and research modes.
  - Advanced flagging and bi-directional LIS connectivity streamline laboratory workflows.
- 

### Clinical Benefits for Hospital Laboratories

- Offers 37 reportable and 34 research parameters for comprehensive diagnostics.
  - Reduces manual intervention with powerful abnormal flagging and autoloader functionality.
  - Delivers reliable, consistent results that support faster clinical decisions and better patient care.
-

## **Emergency Testing Performance**

- STAT sample position prioritizes emergency and critical tests without delay.
  - External compressor and robust hardware ensure uninterrupted operation even in demanding conditions.
  - High throughput and user-friendly interface mean rapid processing during peak hours.
- 

## **Why Hospitals Choose**

- Durable construction and advanced air-handling offer long service life and low maintenance.
  - Automated cleaning systems and disposable tips help prevent carryover and contamination, supporting patient safety.
  - Effortless reagent loading and intuitive touchscreen controls boost user convenience and reduce training time.
- 

## **Transform Your Laboratory Today**

Upgrade your laboratory with the Agappe Mispa HX88 and achieve superior diagnostic capabilities, improved turnaround times, and scalable performance to meet future healthcare demands.

---