

Product Name: Nebulizer New Model

Sr. No.	Description	Amkay
1	Product Name	Nebulizer
2	Set Contain	1) Compressor 2) Adult Mask 3) Child Mask 4) Medicine Bottle Cap with integrated baffle 5) Medicine Bottle Cup 6) Air tube 7) Air Filter
3	Specification of Device	200.49 X 170.08 X 110.94 mm (Length x Width X Height)
4	Cable wire Length	1.40 Meter
5	Material	ABS Plastic (Device) Masks (Plastic)
6	Individual Weight Of Nebulizer	1.556 Kg
7	Technical Data (As Mentioned in User Manual)	Power : 230V AC. 50 Hz Power Consumption : 75 W Medication Capacity : 2ml-6ml Particle Size : 0.4 to 8 um MMAD : 4 um On Time : 30 Minutes Off Time : 40 Minutes Sound Level : 70dB(1 meter) Nebulization rate : 0.2 ml/min or above Compressor Pressure range : 28 PSI Operating Pressure range: 8 to 16 PSI Compressor air output : 10 LPM Operating temperature range: 0°C to 40°C Operating humidity range : 10 to 75% RH Storage humidity range : 90 % RH (Max)

8 Image Of Set



9 Images

